

TECHNICAL SPECIFICATION

Energy adjustment	30-75 speed , continuously adjustment
Step Length	5 speed
Working frequency	230Khz±20Khz
Shots	525 shots, each shot shoot 576 points
Output power	Not more than 580VA
Normal operating conditions	
Ambient Temperature Range	5C-40C
Relatice humidity range	≤80%
The Atmospheric pressure range	86. 0kp-106. 0kpa
Power Supply	AC 110-230V, 50HZ/60HZ

BEFORE&AFTER



Head Office – Sydney
 Suite 1A, Level 2
 802 Pacific Highway
 Gordon NSW 2072
 Tel: (02) 8231 6529

Brisbane Office
 Level 19, 1 Eagle Street
 Brisbane 4000
 Tel: (07) 3535 1202

Melbourne Office
 Level 9, 440 Collins Street
 Melbourne 3000
 Tel: (03) 8618 6934



UNIVERSAL MEDICAL AESTHETICS
 Technology designed for your business

LIPOSUNIC SLIMMING MACHINE-HIFU-II



TECHNICAL FEATURES

The Liposonic system uses high-intensity focused ultrasound delivering uniform energy at a precise depth (13mm) to destroy targeted subcutaneous fat by causing a rapid local temperature rise. This action eliminates targeted adipose tissue just beneath the skin in the abdomen and flanks. Your clients may experience predictable results after a single one-hour treatment. It can penetrate the skin layers and reach the targeted fat tissue without harming the skin or surrounding tissues. Meanwhile the metabolized contents (triglyceride, fatty acid) of cells will excrete out of the body automatically via lymphatic drainage.

The Liposonic treatment uses advanced ultrasound technology to permanently destroy unwanted abdominal fat. Some people may be candidates for the lipohifu treatment of other body areas, including the arms, thighs, and back.

TREATMENT AREA



A Liposonic treatment of the abdomen and flanks takes approximately one hour. After the treatment area is marked, the treatment grid is applied and numbered. The handpiece is placed over the first site and treatment begins. This is repeated until every site has been treated with the appropriate energy.

TREATMENT STEP



1. The treatment head of the Liposonic system is placed on the area to be treated.
2. High-intensity ultrasound energy focused on the targeted fat, destroying fat cells beneath the surface layers of the skin.
3. A continuous lesion is created in the targeted fat and skin layer is not harmed.
4. And where it goes? The treated fat tissue is removed through the body's natural healing and elimination process, leaving you with a thinner, more contoured waistline.

WHAT MAKES IT STAND OUT?

- Equipped with 2 standard heads (13.0mm, 8.0mm)
- Results may be visible after one treatment.
- No damage to the skin.
- No down time, resume your normal activities immediately.
- Dermal collagen stimulation.

